

Statement

To whom it may concern,

In response to the recent updated multiple variants of the novel coronavirus, we, Guangdong Longsee Biomedical Co.,Ltd., as a manufacturer of 2019-nCoV Ag Rapid Detection Kit (Immuno-Chromatography)(LS-C-T-008/009), hereby declare that these products remain effective for the detection of SARS-CoV-2antigen in discovered variants. The variants include but not limited to the novel coronavirus variants such asAlpha (B.1.1.7),Beta (B.1.351/B.1.351.2/B.1.351.3),Gamma(P,1/P.1.1/P,1.2), Delta(B.1.617.2/AY.1/AY.2)Eta(B.1.525),Iota(B.1.526),Kappa(B.1.617.1),Lambda(C.37),Zeta(P.2),Theta(P.3),Epsilon (B.1.427/B.1.429), Omicron(B. 1.1.529/BA.2.12.1/BA.4/BA.5) ,Omicron-Variant EG.5 and EG.5.1,etc. The detection target of our antigen test kits is the nucleocapsid protein (nucleoprotein or protein N) antigens, but these mutations of variants occurred in the spike protein and hardly affect nucleocapsid protein. Therefore, we confirm our products can detect these variants.

Thank you for choosing LONGSEE products, we will continue to track and update the latest viral information including mutations and substitutions for SARS-CoV-2 to upgrade the specificity and sensitivity of our products and strengthen our research and development to provide you with high quality products and services

Guangdong Longsee Biomedical Co., Ltd.

August 31, 2023

